



Product Specification Sheet

TS/18-B

3 inch swab with a flexible polypropylene shaft with cotton tipped swab in 3 inch polypropylene outer tube

Product Code	TS/18-B
Description	3 inch swab with a flexible polypropylene shaft with cotton tipped swab in 3 inch polypropylene outer tube
Intended Use	Ideal for general swabbing where a small shaft is required.
Components	Polypropylene Tube Polystyrene cap cotton tip/polypropylene shaft swab (tip diameter – 5.15 mm) Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	
Pack Size	1000 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.