



Product Specification Sheet

TS/5-6

5ml Amies with charcoal transport medium in tube(black cap). Includes aluminium shaft, cotton mini-tip swab (orange cap) in a peel pouch.

| | |
|-----------------------------|--|
| Product Code | TS/5-6 |
| Description | 5ml Amies with charcoal transport medium in tube(black cap). Includes aluminium shaft, cotton mini-tip swab (orange cap) in a peel pouch. |
| Intended Use | Bacteriological samples only. Collection and transport of clinical specimens |
| Components | Black cap Polypropylene Tube Orange cap cotton/aluminium swab Non-toxic glue Label & film |
| Contents | De-ionised water Amies Charcoal Transport Medium |
| Method of Sterilisation | Sterilised by gamma irradiation at dosage of 25-46kGy |
| Shelf Life | 24 months from date of manufacture |
| Storage Conditions | 4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container |
| Precautions and Limitations | Recovery after 48 hours will be unpredictable and target organisms will be lost. Fatty acid sensitive strains will give decreased recovery or be stressed in the presence of cotton. |
| Pack Size | 100 in a box |
| Other | Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism Recovery |

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 6 Date : 06.06.2016

ProductSpec_TS5-6