



Product Specification Sheet

TS/15-F40

Sterile 500 ml Empty Sampling Container

Product Code	TS/15-F40
Description	Sterile 500 ml Empty Sampling Container
Intended Use	Sampling of potable and swimming pool water samples for microbiological analysis or for general purpose applications.
Components	HDPE Bottle, HDPE Cap, Label.
Contents	Not applicable as supplied empty.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	400 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 3 Date: 24.04.2020

ProductSpec_TS15-F40